


SysTM™ Instrument Cassettes

For efficient cleaning, sterilization and storage of instruments



SciCan
A HIGHER STANDARD



Why use cassettes?

- Virtually eliminates hand scrubbing when used in conjunction with Hydrim® instrument washer
- Efficient, time-saving transfer of instruments from the operator to Hydrim to Statim®
- Keeps instruments organized and neat
- Helps prevent damage to instruments

Features and Benefits of SciCan's Sys™ cassettes

- Unique hole pattern provides excellent washing, rinsing and draining in the Hydrim instrument washer/washer-disinfector
- High grade surgical stainless steel with a passivated finish to maximize durability and minimize discoloration
- Manufacturer's four year warranty on all moving parts and a lifetime guarantee against corrosion



Statim 2000/5000 Cassettes

Select from four sizes, designed to fit the Statim 2000 and 5000.



- SysTM 1 2.5 x 7 x 0.5" (64 x 178 x 13 mm)
- ideal for exam set-ups
 - can hold up to 7 instruments



- SysTM 3 4.5 x 7 x 0.5" (114 x 178 x 13 mm)
- ideal for perio/hygiene
 - capacity of 10 to 12 instruments



- SysTM 2 2.5 x 7 x 1" (64 x 178 x 25 mm)
- perfect size for Statis or other handpieces
 - holds 1 to 2 handpieces, or various other instruments



- SysTM 4B 10.5 x 6.5 x 1" (267 x 165 x 25 mm)
- designed for larger operative/surgical set-ups
 - capacity of 14 instruments plus an accessory area



- Modular design maximizes space utilization and optimizes the efficiency of your Statim sterilizer



- Silicone instrument retention strips are designed to minimize contact area with instruments, which promotes more effective cleaning
- Choice of 8 colours of silicone strips to customize cassettes for your needs (all cassettes come standard with green silicone colour coding)

Statim 7000 Cassettes

Select from five available sizes, designed to fit the Statim 7000 in a modular fashion, and ensure effective use of space. Double hinged for easy lid management at chairside. Scalloped instrument supports keep instruments in place.



Full-size 8 x 11 x 1" (203 x 280 x 25 mm)

- designed for larger surgical set-ups
- capacity of 20 instruments, or 14 instruments and accessory area
- side-mounted latches for easy opening



Half-size 5.5 x 8 x 1" (140 x 203 x 25 mm)

- ideal for hygiene, perio, prophy
- capacity of 10 instruments



Half-size, Double Stack 5.5 x 8 x 1.5" (140 x 203 x 38 mm)

- stores twice the instruments as the standard half size, in the same footprint
- capacity of 20 instruments



Quarter-size 2.75 x 8 x 1" (69 x 203 x 25 mm)

- perfect for exam
- capacity of 5 instruments



Quarter-size, Double Stack 2.75 x 8 x 1.5" (69 x 203 x 38 mm)

- stores twice the instruments as the standard quarter size, in the same footprint
- capacity of 10 instruments

Cassette Part-number	Cassette Description	Instrument Capacity	Max. Cassette Load**		
			Statim* 7000	Hydrim* L110w/M2	Hydrim* C51w/wd
Sys-TM-8 x 11	Full-size	20	2	10	4
Sys-TM-8 x 11A	Full-size with accessory area	14	2	10	4
Sys-TM-8 x 5.5	Half-size	10	4	20	8
Sys-TM-8 x 5.5DS	Half-size, double stack	20	2	20	8
Sys-TM-8 x 2.75	Quarter-size, double stack	5	8	40	16
Sys-TM- 8 x 2.75DS	Quarter-size, double stack	10	4	40	16

Cassette Part-number	Cassette Description	Instrument Capacity	Max. Cassette Load**			
			Statim* 2000/2000S	Statim* 5000/5000S	Hydrim* L110w/M2	Hydrim* C51w/wd
Sys-TM1	Exam	7	6	10	40	16
Sys-TM2	Handpiece/ surgery	Varies	4	8	40	16
Sys-TM3	Perio/hygiene	10-12	4	8	20	12
Sys-TM4B	Operative/ surgical	14	1	2	10	4

Colour Coding Part-number	Description
VSTSB-PSG	1' (30 cm) Silicone Cassette Rail - Blue
VSTSG-PSG	1' (30 cm) Silicone Cassette Rail - Green
VSTSK-PSG	1' (30 cm) Silicone Cassette Rail - Black
VSTSL-PSG	1' (30 cm) Silicone Cassette Rail - Light Blue
VSTSP-PSG	1' (30 cm) Silicone Cassette Rail - Pink
VSTSR-PSG	1' (30 cm) Silicone Cassette Rail - Red
VSTSU-PSG	1' (30 cm) Silicone Cassette Rail - Purple
VSTSY-PSG	1' (30 cm) Silicone Cassette Rail - Yellow



www.scican.com

SciCan Ltd., 1440 Don Mills Road, Toronto Ontario, Canada M3B 3P9
Telephone: 1.800.667.7733 Fax: 416.445.2727

SciCan Inc., 701 Technology Drive, Canonsburg, PA 15317
Telephone: 1.800.572.1211 Fax: 724.820.147

EU Representative: SciCan GmbH, Kurzes Geländ 10, 86156 Augsburg, Germany
Telephone: +49 (0)821.56 74 560
Fax: +49 (0)821.56 74 56 99

SciCan Medtech, Alpenstrasse 16, 6300 Zug, Switzerland
Telephone: +41.41.727.7027 Fax: +41.41.727.7029

Hydrim and Statim are registered trademarks of SciCan Ltd. SysTM is a trademark of SciCan Ltd.
* Not all Statim and Hydrim models are available in all markets.

SD-294-E 10M 9-08

SciCan

A HIGHER STANDARD

**Based on empty cassette weight. Operator must adhere to weight restrictions for full loads as listed in the sterilizer operator manual.